ENGLISH

E

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION

IMPORTANT

IThis booklet is designed to assist in using the D E N A LI®, DENALI® Deformity, MESA® and RANGE® Spinal Systems. It is not a reference for surgical techniques.

CAUTION: Federal law (USA) restricts this device to sale and use by, or on the order of, a physician.

INDICATIONS

The RANGE Spinal System is comprised of the DENALI, DENALI Deformity and MESA Spinal Systems all of which are cleared for the following indications:

Non-cervical spinal fixation devices intended for posterior, non-pedicle stabilization as an adjunct to fusion for the following indications: degenerative disc disease (DDD) (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e. fracture or dislocation); spinal stenosis; curvatures (i.e. scoliosis, kyphosis; and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion.

Non-cervical pedicle screw fixation devices for stabilization as an adjunct to fusion for the following indications: trauma (i.e. fracture or dislocation); spinal stenosis; curvatures (i.e. scoliosis, kyphosis; and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion. It is also indicated for the treatment of severe spondylolisthesis (grades 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone graft having implants attached to the lumbar and sacral spine (L3 to sacrum) with removal of the implants after the attainment of a solid fusion.

MATERIALS

All implant components are manufactured from Titanium alloy, CP Titanium and Cobalt Chrome per ASTM and ISO standards.

CLEANING/ REPROCESSING OF K2M SURGICAL INSTRUMENTS

K2M surgical instruments are supplied non-sterile. While it is recommended that the following steps are included in a decontamination/ reprocessing protocol the end-user bears the ultimate responsibility for the cleanliness of the device. These instructions are not intended for K2M implants or disposable surgical instruments.

Presoak the instruments with an enzymatic solution for a minimum of 5 minutes. Following the presoak, the instruments should be wiped or scrubbed using a brush, cloth or sponge that does not mar the surface of the instrument. Remove soil from cannulated parts with a nylon bristle brush or appropriately sized guide wire. Rinse parts under water for one minute. Repeat the process until no visible debris remains. Clean K2M surgical instruments with an appropriate brush, cloth or sponge and low foaming, pH neutral detergent solution. The use of abrasive compounds or excessively acidic or alkaline solutions may cause damage to the instruments and should be avoided. Rinse parts under warm or hot flowing water for a minimum of 1 minute including direct contact with all surfaces for at least 10 seconds. Repeat rinsing step using distilled, reverse osmosis or deionized water. Automatic cleaning may be used in addition to manual cleaning. Do not ultrasonically clean tools requiring limiting handles.

STERILIZATION

Unless specifically labeled sterile, the implants and instruments are supplied NONSTERILE and MUST be sterilized prior to use. Recommended sterilization methods include steam autoclaving after removal of all protective packaging and labeling. The following steam autoclave cycles are recommended, however sterilization should be in accordance with the sterilizer manufacturer’s instructions and the institution’s procedures for assuring sterility.

AUTOClave CYCLES:

Prevacuum

Temperature: 270°F-275°F (132°C-135°C) Time: 20 minutes

Use caution during sterilization and storage. Do not allow contact with metal or other hard objects that could damage the finish or prevent proper use. (See Preoperative Warnings and Precautions).

NOTE: Instruments that may have been exposed to Creutzfeldt-Jakob disease (CJD) should be treated according to the hospital’s prion decontamination protocol. K2M recommends contacting the Centers for Disease Control and the World Health Organization for the most recent information on CJD transmission and deactivation.

INSTRUCTIONS FOR USE [ⓘ]

(For complete instructions refer to the appropriate surgical technique provided by your local K2M sales representative.)

The patient is placed in the position desired by the surgeon to allow a standard approach to the appropriate portion of the patient’s anatomy for the procedure.

Following a standard approach to the patient’s spine the appropriate implants are used for either screw fixation or hook fixation.

NOTE: Excessive reverse bending of Titanium Rods can cause metal stressing resulting in a lower fatigue life for the rod.

Screw Fixation

For screw fixation use the probe, reamer and tap to prepare the screw site. Select the proper size screw. Insert with the screw inserter.

NOTE: Taps are sized smaller than the actual diameter of the screw to allow better screw fixation.

Cut the rod to the proper length. Bend the rod as needed due to anatomical variations. Drop the rod into the top of the screw.

Complete final locking of the assembly as described in the appropriate surgical technique.

Once the surgeon is satisfied the device has been properly implanted, the surgical site is closed in the usual manner.

For Hook Use

Prepare the hook site with the appropriate lamina finders to provide a good fit with the hooks. Assemble the proper hook onto the hook holder. Insert hook using the hook holder. Remove the hook holder.

Cut the rod to the proper length. Bend the rod as needed due to anatomical variations. Drop the rod into the top of the hook.

Insert the set screw. Use distraction or compression provided by the distractor or compressor as required. Thread the set screw completely down across the rod the assembly.

Once the surgeon is satisfied the device has been properly implanted, the surgical site is closed in the usual manner.

CONTRAINDICATIONS

1. K2M spinal systems are contraindicated in the presence of infection, pregnancy, metabolic disorders of calcified tissues, grossly distorted anatomy, inadequate tissue coverage, drug/alcohol abuse, mental illness, general neurological conditions, immunosuppressive disorders, patients with known sensitivity to materials in the device, obesity, patients who are unwilling to restrict activities or follow medical advice, and any condition where the implants interfere with anatomical structures or precludes the benefit of spinal surgery.

2. Biological factors such as smoking, use of nonsteroidal anti-inflammatory agents, the use of anticoagulants, etc. all have a negative effect on bony union. Contraindications may be relative or absolute and must be carefully weighed against the patient’s entire evaluation.

3. This device is not intended for use except as indicated.

POTENTIAL ADVERSE EVENTS

1. Potential adverse events include, but are not limited to pseudoarthrosis; loosening, bending, cracking or fracture of components, or loss of fixation in the bony union with possible neurologic damage, usually attributable to pseudoarthrosis, insufficient bone stock, excessive activity or lifting, or one or more of the factors listed in Contraindications; or Warnings and Precautions; infections possibly requiring removal of devices; palpable components, painful fusions, and/or pressure necrosis; and allergies; and other reactions to device materials which, although infrequent, should be considered, tested for (if appropriate), and ruled out pre-operatively.

2. Potential risks also include those associated with any spinal surgery resulting in neurological, cardiovascular, respiratory, gastrointestinal or reproductive compromise, or death.

WARNINGS AND PRECAUTIONS

1. The K2M Pedicle Screw and Hook Systems are intended for use for the indications listed. Safety and effectiveness of the implants have not been established for other applications. The implants are for single use only and are not designed to be combined with devices from other manufacturers.

2. For optimum results careful preoperative diagnosis and planning, meticulous surgical technique and extended postoperative care by experienced spinal surgeons are essential. Prior to use the surgeon should be specifically trained in the use of this spinal system and the associated instrumentation to facilitate correct selection and placement of the implants. The size and shape of bones and soft

tissue place limitations on the size and strength of the implants and proper selection will reduce the risk of neurological injury during implantation as well as metal fatigue leading to bending or breakage of the device.

3. Patient selection and compliance is extremely important. Spinal implant surgery on patients with conditions listed under Contraindications may not be candidates for this procedure. The patient must be made aware of the limitations of the implant and that physical activity and load bearing have been implicated in premature loosening, bending or fracture of internal fixation devices. The patient should understand that a metallic implant is not as strong as a normal, healthy bone and will fracture under normal load bearing in the absence of complete bony healing. An active, debilitated/ uncooperative patient who cannot properly restrict activities may be a particular risk during postoperative rehabilitation.

4. Potential risks identified with the use of this device system which may require additional surgery include device component failure, loss of fixation, non-union, fracture of the vertebra, and neurological, vascular or visceral injury.

5. Cutting, bending, or scratching the surface of metal components can significantly reduce the strength and fatigue resistance of the implant system and should be avoided where possible. These, in turn may cause cracks and/or internal stresses that are not obvious to the eye and may lead to fracture of the components. Especially avoid sharp or reverse bends and notches.

6. Special protection of implants and instruments during storage is recommended when exposed to corrosive environments such as moisture, salt, air, etc.

7. Implanting metals and alloys in the human body subjects them to a constantly changing environment of salts, acids and alkalis which can cause corrosion. Putting dissimilar metals in contact with each other can accelerate the corrosion process which in turn may enhance fatigue fractures of implants. Thus every effort should be made to use compatible metals and alloys. Fretting or wear at the interface between components of a device may also accelerate the corrosion process and may lead to the generation of wear debris which has been associated with localized inflammatory response.

8. The K2M spinal implants are intended to provide temporary stabilization. If an implant remains implanted after complete healing it can actually increase the risk of refracture in an active individual. The surgeon should weigh the risks versus the benefits when deciding whether to remove the implant.

PREOPERATIVE

1. Patient conditions and/or predispositions such as those previously addressed in Contraindications and Warnings and Precautions should be avoided.

2. Preoperative testing (simple bend and where necessary, stretch testing) should identify degree of correction possible without neurological damage and levels to be spanned using techniques similar to other spinal fusion procedures.

3. Use care in handling and storage of the implants. Prior to surgery components should be inspected for any evidence of damage or corrosion.

4. An adequate inventory of implant sizes should be available at the time of the surgery.

5. All components should be cleaned and sterilized before use.

6. Before the initial experience we recommend that the surgeon critically review all available information and consult with other surgeons having experience with the device.

OPERATIVE

1. The primary goal of this surgery is to arthrodesis selected vertebrae. Adequate exposure, bony preparation and grafting are essential to achieving this result.

2. Rods may be prebent to the degree of correction determined by preoperative testing however reverse bends should be avoided.

3. The use of two rods and crosslinking the rods will provide a more rigid construct.

4. The placement of screws should be checked radiographically prior to assembly of the rod construct.

5. Care should be taken when positioning the implants to avoid neurological damage.

6. Use of bone cement will make removal of the implants difficult and should be avoided.

POSTOPERATIVE

1. Adequately instruct the patient. Postoperative care and the patient’s ability and willingness to follow instructions are two of the most important aspects of successful healing.

2. Internal fixation devices are load sharing devices which maintain alignment until healing occurs. If healing is delayed or does not occur the implant could eventually break, bend or loosen. Loads produced by load bearing and activity levels will impact the longevity of the implant.

3. Metallic implants can loosen, fracture, corrode, migrate, cause pain, or stress shield bone even after a bone has healed. If an implant remains implanted after complete healing, it can actually increase the risk of refracture in an active individual. The surgeon should weigh the risks versus benefits when deciding whether to remove the implant. Implant removal should be followed by adequate postoperative management to avoid refracture.

4. Periodic X-rays for at least the first year postoperatively are recommended for close comparison with postoperative conditions to detect any evidence of changes in position, nonunion, loosening, and bending or cracking of components. With evidence of these conditions, patients should be closely observed, the possibilities of further deterioration evaluated, and the benefits of reduced activity and/or early revision considered.

5. Surgical implants must never be reused. An explanted metal implant should never be reimplanted. Even though the device appears undamaged, it may have small imperfections and internal stress patterns which may lead to early breakage.

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